








## Dräger Savina® 300 Disassembly/Assembly



### Disassembly/Assembly Flow Sensor

<p><b>1</b> Open the flap by lifting the lower edge upwards.</p>  <p style="text-align: right;"><b>3</b> Close the flap.</p>	<p><b>2</b> Push the flow sensor as far as it will go to the left.</p>  <p style="text-align: right;"><b>2</b> Push the flow sensor to the right into the socket.</p>	<p><b>3</b> Remove the flow sensor from the socket.</p>  <p style="text-align: right;"><b>1</b> Fit the flow sensor into the socket.</p>
--	---	--

### Disassembly/Assembly Expiratory Valve

<p><b>1</b> Turn the locking ring as far as possible to the left.</p>  <p style="text-align: right;"><b>5</b> Turn the locking ring as far as possible to the right until it clicks into place.</p>	<p><b>2</b> Remove the expiratory valve from the fitting.</p>  <p style="text-align: right;"><b>4</b> Turn the locking ring as far as possible to the left, push the expiratory valve into the fitting.</p>	<p><b>3</b> Pull the sensor sleeve from the expiratory valve.</p>  <p style="text-align: right;"><b>3</b> Fit the flow sensor sleeve on the expiratory valve.</p>
<p><b>4</b> Remove the diaphragm.</p>  <p style="text-align: right;"><b>2</b> Fit the diaphragm onto the edge of the expiratory valve housing.</p>	<p><b>5</b> Remove and empty the water trap container.</p>  <p style="text-align: right;"><b>1</b> Fit the water trap container.</p>	

This guideline does not replace or amend the instructions for use.  
The instructions for use and the hospital guidelines and policies always have to be observed.

## Dräger Savina® 300

### Flow Sensor: Reprocessing instructions Manual cleaning and disinfection



**1** Prepare cleaning solution.



Immerse flow sensor bubble-free in the cleaning solution for the agent specific exposure time.

**2** At the beginning and at the end vigorously swirl the flow sensor in the cleaning solution at least 3 times.



Ensure that the cleaning solution flushes all surfaces and interior spaces.

**3** Rinse the flow sensor in a water (minimum drinking-water quality) and wash off all cleaning agent residues.



**4** Prepare the disinfectant solution.




Repeatedly swirl the flow sensor in the disinfection solution for the agent specific exposure time. Ensure that the solution flushes all surfaces and interior spaces.

**5** Sufficiently rinse the flow sensor in water (minimum drinking-water quality) until no disinfectant residues can be recognized.



Inspect flow sensor for visible soiling and damage.

**6** Thoroughly shake off remaining water.



Allow the flow sensor to air dry for at least 30 minutes.

In case of a combined cleaning-disinfection agent please also follow the guideline as described (two times immersion).


## Dräger Savina® 300

Reprocessing instructions

Manual cleaning and disinfection



### Reprocessing instructions overview:

	Manual cleaning and disinfection	Machine cleaning and disinfection	Sterilization
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✗	✓
	✓	✗	✗

#### Machine cleaning and disinfection:

Securely position the product in the basket.

Ensure the following:

- All surfaces and interior spaces can be flushed completely.
- The water can drain off freely.

#### Sterilization:

Check the product for visible damage and replace if necessary.

For these semi critical products sterilization is possible but not mandatory.

**For further information please refer to Instruction for Use.**

This guideline does not replace or amend the instructions for use.

The instructions for use and the hospital guidelines and policies always have to be observed.

#### Manufacturer:

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany